

FAME Study

FFR vs Angiography for Multivessel Evaluation



ST. JUDE MEDICAL


MORE CONTROL. LESS RISK.

Does measuring FFR really make a difference?

FAME is a large, randomized, controlled, multicenter trial comparing stent treatment guided by FFR versus angiographic guidance alone

The FAME study was designed to reflect the daily practice of performing PCI in patients with multivessel disease

Study Design



Key Inclusion and Exclusion criteria

Inclusion criteria:

- Patients with multivessel disease
- At least 2 stenoses greater than or equal to 50% in 2 or 3 major epicardial coronary artery vessels, amenable to stenting

Exclusion criteria:

- Left main disease or previous bypass surgery
- ST-elevation MI with CK greater than 1000 U/l within last 5 days
- Extremely tortuous or calcified coronary arteries

Note: patients with previous PCI were not excluded



Endpoints

- **Primary Endpoint**

Composite of death, myocardial infarction, or repeat revascularization (“MACE”) at 1-year

Secondary Endpoints

- individual components of MACE at 1-year
- functional class
- use of anti-anginal drugs
- health-related quality of life (EuroQOL-5D)
- procedure time
- amount of contrast agent used during procedure
- cost of the procedure

Participating Sites

Fourteen European centers and six United States centers enrolled more than 1,000 patients in the FAME Study

September 25, 2007


Total Enrollment Concluded

1005

Inclusion Completed



Study Enrollment and Randomization



Baseline Characteristics

Characteristic	Angiography Group (N= 496)	FFR Group (N= 509)	P Value†
Demographic			
Age — yr	64.2±10.2	64.6±10.3	0.47
Sex — no. (%)			0.30
Male	360 (72.6)	384 (75.4)	
Female	136 (27.4)	125 (24.6)	
Clinical			
Angina classification — no. (%) †			0.13
I	115 (23.2)	132 (25.9)	
II	165 (33.3)	170 (33.4)	
III	118 (23.8)	132 (25.9)	
IV	98 (19.8)	75 (14.7)	
Previous myocardial infarction — no. (%)	180 (36.3)	187 (36.7)	0.84
Previous PCI — no. (%)	129 (26.0)	146 (28.7)	0.34
Diabetes — no. (%)	125 (25.2)	123 (24.2)	0.65
Hypertension — no. (%)	327 (65.9)	312 (61.3)	0.10
Hypercholesterolemia — no. (%)	362 (73.0)	366 (71.9)	0.62
Family history — no. (%)	190 (38.3)	205 (40.3)	0.49
Current smoker — no. (%)	156 (31.5)	138 (27.1)	0.12
Unstable angina — no. (%)			
With dynamic ECG changes	91 (18.3)	73 (14.3)	0.09
Without dynamic ECG changes	87 (17.5)	77 (15.1)	0.29
Left ventricular ejection fraction — %	57.1±12.0	57.2±11.0	0.92
Medication			
Beta-blocker — no. (%)	377 (76.0)	395 (77.6)	0.55
Calcium antagonist — no. (%)	96 (19.4)	121 (23.8)	0.09
Nitrate — no. (%)	179 (36.1)	167 (32.8)	0.27
ACE inhibitor or ARB — no. (%)	255 (51.4)	267 (52.5)	0.74
Statin — no. (%)	397 (80.0)	417 (81.9)	0.45
Aspirin — no. (%)	454 (91.5)	465 (91.4)	0.92
Clopidogrel — no. (%)	292 (58.9)	310 (60.9)	0.51

Similar characteristics in the two groups



Baseline Characteristics (cont'd)

Characteristic	Angiography Group (N=496)	FFR Group (N=509)	P Value†
Angiographic Findings			
Indicated lesions per patient — no.‡	2.7±0.9	2.8±1.0	0.34
Extent of occlusion — no. of lesions/total no. (%)			
50–70% narrowing	550/1350 (40.7)	624/1414 (44.1)	
71–90% narrowing	553/1350 (41.0)	530/1414 (37.5)	
91–99% narrowing	207/1350 (15.3)	202/1414 (14.3)	
Total occlusion	40/1350 (3.0)	58/1414 (4.1)	
Patients with total occlusion — no. (%)	37 (7.5)	54 (10.6)	
Quantitative coronary analysis			
Extent of stenosis — %	61.2±16.6	60.4±17.6	0.24
Minimal luminal diameter — mm	1.0±0.4	1.0±0.5	0.35
Reference diameter — mm	2.5±0.6	2.5±0.7	0.81
Lesion length — mm	12.6±6.9	12.5±6.5	0.42
SYNTAX score¶	14.5±8.8	14.5±8.6	0.95
EQ-5D score	64.7±19.2	66.5±18.3	0.24

Results

Variable	Angiography Group (N=496)	FFR Group (N=509)	P Value†
Procedure time — min‡	70±44	71±43	0.51
Volume of contrast agent used — ml	302±127	272±133	<0.001
Drug-eluting stents			
No. of stents per patient			
Mean	2.7±1.2	1.9±1.3	<0.001
Median (interquartile range)	3 (2–3)	2 (1–3)	
Total length per patient — mm	51.9±24.6	37.9±27.8	<0.001
Average diameter per patient — mm	2.96±0.33	2.92±0.36	0.13
Total no. of stents			
Zotarolimus-eluting — no. (%)	603 (44.4)	403 (41.1)	
Sirolimus-eluting — no. (%)	273 (20.1)	202 (20.6)	
Paclitaxel-eluting — no. (%)	414 (30.5)	316 (32.2)	
Other — no. (%)	69 (5.1)	59 (6.0)	
Lesions in which stents successfully placed — no./total no. (%)§	1237/1350 (91.6)	819/874 (93.7)	
FFR-guided strategy			
Lesions successfully measured for FFR — no./total no. (%)¶	NA	1329/1414 (94.0)	
FFR	NA	0.71±0.18	
Ischemic lesions	NA	0.60±0.14	
Nonischemic lesions	NA	0.88±0.05	
Lesions with FFR ≤0.80 — no./total no. (%)	NA	874/1387 (63.0)	
Lesions with FFR >0.80 — no./total no. (%)	NA	513/1387 (37.0)	
Cost of materials — \$	6,007±2,819	5,332±3,261	<0.001
Hospital stay at baseline admission — days	3.7±3.5	3.4±3.3	0.05

**Significant difference
between the two
groups**



Endpoints

End Point	Angiography Group (N= 496)	FFR Group (N= 509)	P Value†	Relative Risk with FFR Guidance (95%CI)
Events at 1 year				
Composite of death, myocardial infarction, and repeat vascularization — no. (%)‡	91 (18.3)	67 (13.2)	0.02	0.72 (0.54–0.96)
Death — no. (%)	15 (3.0)	9 (1.8)	0.19	0.58 (0.26–1.32)
Myocardial infarction — no. (%)	43 (8.7)	29 (5.7)	0.07	0.66 (0.42–1.04)
Repeat vascularization — no. (%)	47 (9.5)	33 (6.5)	0.08	0.68 (0.45–1.05)
Death or myocardial infarction — no. (%)	55 (11.1)	37 (7.3)	0.04	0.66 (0.44–0.98)
Total events — no.	113	76		
Events per patient — no.	0.23±0.53	0.15±0.41	0.02	
Functional status at 1 year				
Patients without event and free from angina — no./total no. (%)	326/482 (67.6)	360/493 (73.0)	0.07	
Patients free from angina — no./total no. (%)	374/480 (77.9)	399/491 (81.3)	0.20	
Antianginal medications — no.§	1.23±0.74	1.20±0.76	0.48	
Score on EQ-5D visual-analogue scale¶	73.7±16.0	74.5±15.7	0.65	



Results One-year Follow-up

FAME:

- proof that measuring FFR during the stenting procedure really does make a difference...

Results One-year Follow-up

...there is a reduction in MACE: 28%

less risk of dying, having a heart attack or having to come back for more stents or a bypass operation

...there is a reduction in death or MI: 34%


less risk of dying or having a heart attack

Summary:

**Better outcomes... AND it significantly SAVES on costs
AND it doesn't prolong the procedure.**



Sub-analysis – Angiographic vs Functional Severity



Tonino P. A. L et al; *J. Am. Coll. Cardiol.* 2010;55;2816-2821



ST. JUDE MEDICAL™

MORE CONTROL. LESS RISK.


Results 2-year Follow-up

Two-year Results, No. (%)	ANGIO-Group N=496	FFR-Group N=509	P-Value
Death, MI and Repeat Vascularization	110 (22.2)	90 (17.7)	0.07
Death or Myocardial Infarction	63 (12.7)	43 (8.4)	0.03
Death	19 (3.8)	13 (2.6)	0.25
Myocardial Infarction	48 (9.7)	31 (6.1)	0.03
Freedom from Angina	(75.8)	(79.9)	0.14
Follow-up	(92.7)	(94.5)	0.34

Significant difference in MI and MI/death between the two groups

Tonino P. A. L *et al*; *J. Am. Coll. Cardiol.* 2010;55;2816-2821

2-year Survival from Death and MI




Tonino P. A. L *et al*; *J. Am. Coll. Cardiol.* 2010;55;2816-2821

Results Two-year Follow-up

Two-year Results, No. (%)	ANGIO-Group N=496	FFR-Group N=509	P-Value
Death, MI and Repeat Vascularization	110 (22.2)	90 (17.7)	0.07
Death or Myocardial Infarction	63 (12.7)	43 (8.4)	0.03
Death	19 (3.8)	13 (2.6)	0.25
Myocardial Infarction	48 (9.7)	31 (6.1)	0.03
Freedom from Angina	(75.8)	(79.9)	0.14
Follow-up	(92.7)	(94.5)	0.34


Significant difference in MI and MI/death between the two groups

2-year Survival from Death and MI



Tonino P. A. L *et al*; *J. Am. Coll. Cardiol.* 2010;55;2816-2821

Outcome of Deferred Lesions




William F. Fearon et al, on behalf of the FAME Study Investigators.
Presented at TCT “==)



ST. JUDE MEDICAL™

MORE CONTROL. LESS RISK.

Outcome of Deferred Lesions



William F. Fearon et al, on behalf of the FAME Study Investigators.
Presented at TCT “==)



2-year follow-up: Better Outcomes at Lower Costs



Bootstrap simulation indicated that the FFR-guided strategy was cost-saving in 99.8% and cost-effective in all 1,000 scenarios.

Tonino P. A. L et al; *J. Am. Coll. Cardiol.* 2010;55;2816-2821

Measuring FFR using St Jude Medical's PressureWire improves patient outcomes up to 2-years post procedure and reduces procedural *and* healthcare costs without prolonging procedure time.

Results Two-year Follow-up

- At 2-years a significant decrease in the rate of MI in the FFR-guided arm is demonstrated.
- There continues to be a significant decrease in death and MI favoring the FFR-guided approach.
- There is a strong trend toward a lower rate of death, MI or the need for repeat revascularization in the FFR-guided arm.
- There is no sign suggesting that deferred lesions are likely to be responsible for late myocardial infarctions or to progress and require repeat revascularizations.



Rx Only

Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Product referenced is approved for CE Mark.

PressureWire is designed, developed and manufactured by St. Jude Medical Systems AB. PressureWire, RADI, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.

©2011 St. Jude Medical, Inc. All rights reserved.