

RadiAnalyzer™ Xpress

Technical Specifications



ST. JUDE MEDICAL

MORE CONTROL. LESS RISK.

Technical Specifications

PRESSURE	
Range	-30 to 300 mmHg
Accuracy	± 1 mmHg plus $\pm 1\%$ of reading (-30 to 50 mmHg) $\pm 3\%$ of reading (50 to 300 mmHg)
Resolution	0.1 mmHg
Frequency response	0-25 Hz (standard configuration) 0-100 Hz (with PhysioMon Advanced Analysis option)

Technical Specifications

TEMPERATURE (OPTIONAL)	
Range	15-42°C
Accuracy	0.05° or 10% ΔT whichever greatest
Frequency response	0-7 Hz (temperature mode) 0-25 Hz (CFR mode)

Technical Specifications

ELECTRICAL	
Operating voltage	100-240 VAC, 50-60 Hz

DIMENSIONS AND WEIGHT	
Height x depth x width	29 x 12 x 31.5 cm (11.4 x 4.7 x 12.4 in)
Weight	4.4 kg (9.7 lbs) with transformer 3.7 kg (8.1 lbs) without transformer

Technical Specifications

ENVIRONMENT	
Operating temperature	15-30°C (59-86°F)
Operating relative humidity	30-75%

Rx Only

Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Product referenced is approved for CE Mark.

RadiAnalyzer Xpress, is designed, developed and manufactured by St. Jude Medical Systems AB. RadiAnalyzer Xpress, RADI, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved